

Endozime® AW Plus

Item No. 34514

Multi-Tiered Enzymatic Detergent

Specifically Designed for Washer Sterilizers, Ultrasonics and Manual Cleaning.

A unique synergistic formulation of protease, amylase, lipase, carbohydrase and proprietary enzymes that remove fat, blood, protein, starches and carbohydrates from all surgical instruments and scopes.

The ultimate enzymatic detergent designed as a total cleaning product for use in all washers, washer sterilizers, ultrasonics and manual cleaning. Low-sudsing and neutral pH Endozime® AW Plus eliminates the need for all other cleaners or detergents. Developed for cleaning the most difficult instruments (i.e., orthopaedic, laparoscopic), Endozime® AW Plus has proven to be safe for use on all instruments, even the most delicate (i.e., ophthalmic, microsurgical). At 1/2 ounce per gallon of water, Endozime® AW Plus is the most powerful enzymatic cleaner and a complete detergent with buffers and surfactants. One solution for all your cleaning problems. Clinically tested to be free rinsing and 100% biodegradable.



Features

- Proprietary synergistic mixture of enzymes and non-ionic detergents
- Protease, amylase, lipase, carbohydrase and proprietary enzymes
- Instantly active - Constantly active
- Starts cleaning on contact
- Easy and safe to use
- Safe for all instruments and scopes - Extends usable lifetime
- Neutral pH (non-abrasive)
- Non-foaming - Non-sudsing
- 100% Biodegradable
- Free rinsing
- Will not harm any metals, plastics, rubber or corrugated tubing

Directions for Use

Surgical Instruments and Scopes

1. Add 1/2 ounce of Endozime® AW Plus to 1 gallon (4ml/liter) warm water.
2. Submerge instruments and scopes to be cleaned.
For scopes: suction or flush through channels before soaking.
3. Soak for 2 minutes to remove all organic matter.
4. Rinse thoroughly with tap, distilled or sterile water.

Works in all temperatures.

Discard used solution.

Cautions: Avoid prolonged contact with skin. Do not swallow. Keep away from children.

Automatic Cleaning Devices:

1. Pour directly into detergent reservoir full strength, or
2. Dilute 1/2 ounce per 1 gallon (4ml/liter) of water.