

Prestige Medical - integrated decontamination solutions

Prestige Medical Ltd is an ISO 9001-2000/EN13485 quality accredited company. The Classic is designed to the latest international standards for safety and performance and is fully CE marked to the Medical Devices Directive 93/42/EEC.



Contact us for advice on infection control



Prestige Medical Limited
East House, Duttons Way,
Shadsworth Business Park,
Blackburn BB1 2QR

Tel : +44 (0) 1254 682 622
Fax : +44 (0) 1254 682 606

www.prestigemedical.co.uk
sales@prestigemedical.co.uk

[@PrestigeMedica1](https://twitter.com/PrestigeMedica1)

Registered in England
Reg No. 2826793
Vat No. 633707835

Registered Office:
c/o United Cast Bar (UK) Limited
Spital Lane
Chesterfield
Derbyshire
S41 0EX, England



Prestige Medical is a registered trademark

Due to a continuous process of development, Prestige Medical reserves the right to amend the specification of products and services without notice.

V5 Nov 2016

Classic Portable Autoclave





Prestige Medical Classic – with more than 1 million sold in 80 countries worldwide, you can rely on our reputation for quality and performance.



Technical Specification

Classic Portable Autoclave									
Model Number	210001	210002	210004	210006	210011	210047	210048	210052	210134
Max width - overall	340mm	340mm	340mm	340mm	340mm	340mm	340mm	340mm	340mm
Max height - overall	335mm	420mm	420mm	335mm	335mm	335mm	420mm	335mm	335mm
Max net weight	5.6kg	6.1kg	6.5kg	7.0kg	5.95kg	5.8kg	6.44kg	7.94kg	7.1kg
Chamber diameter	210mm	210mm	210mm	210mm	210mm	210mm	210mm	210mm	210mm
Chamber height	236mm	328mm	328mm	236mm	236mm	236mm	328mm	236mm	236mm
Chamber capacity	9 litres	12 litres	12 litres	9 litres	9 litres	9 litres	12 litres	9 litres	9 litres
Max instrument length	228mm	290mm	290mm	228mm	228mm	228mm	290mm	228mm	228mm
Max load weight	3kg	4kg	4kg	3kg	3kg	3kg	4kg	3kg	3kg
Sterilizing temperature	126°C	126°C	126°C	121°C	126°C	121°C	121°C	134°C	134°C
Min sterilizing time	11 minutes	11 minutes	11 minutes	18 minutes	11 minutes	28 minutes	28 minutes	3 minutes	3 minutes
Cycle time - normal	22 minutes	22 minutes	22 minutes	28 minutes	22 minutes	70 minutes	70 minutes	20 minutes	20 minutes
Operating pressure	1.4 bar	1.4 bar	1.4 bar	1.05 bar	1.4 bar	1.05 bar	1.05 bar	2.2 bar	2.2 bar
Voltage	230v	230v	230v	120v	230v	230v	230v	230v	230v
Wattage	1250w	1250w	1250w	1250w	1250w	1250w	1250w	1250w	1250w
Frequency	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz	50-60Hz
Printer								✓	
Temperature & Pressure Gauges			✓						✓
Furniture	219708	219720	219720	219708	219708	219708	219720	219708	219708
					219706				
					219294				

All Classic units are supplied complete with electrical cord set and instruction manual as appropriate.

Reliable and trusted worldwide

Designed to give you the most cost effective sterilization solution with a wide choice of models and accessories that are guaranteed to meet your needs.

Steam sterilization, or autoclaving, is the preferred method of eliminating the risk of cross infection through contaminated instruments. It is the only method recommended throughout the world by healthcare professionals including the World Health Organisation (WHO) and the Centre for Disease Control (CDC).