

IV
Solutions

Enteral
Care

Respiratory
Care

Medical/
Surgical

Urology

The **Dale Bendable ArmBoard** positions and stabilizes the wrist to prevent unwanted movement when using arterial and peripheral lines in the critical care units, operating room, recovery room, and patient care units.

Dale[®]
Always Better

Bendable ArmBoard

THE ARMBEARD CAN BE CUSTOM-SHAPED TO ANY DESIRED POSITION.

The ArmBoard can be shaped to achieve more prominent radial artery access. A slight bend of the ArmBoard will provide the patient with a more natural hand/wrist position. The small and x-small sizes can be shaped to accommodate an infant's wrist or heel. The ArmBoards are covered in soft, comfortable material, and fastened to the wrist/heel with adjustable straps with hook and loop closure which allows the clinician easy access to the catheter site.



- Bendable inner core
- Comfortably maintains wrist in hyperextended position
- Secures the wrist
- Improves skin care
- Improves patient comfort
- Trouble-free monitoring
- Not made with natural rubber latex.



MR Conditional, See Conditions For Use by visiting our website for detailed information on scanning safely with this device.

ORDERING INFORMATION

Product Number	Size	Quantity
650	Large, 23cm x 9cm (9" x 3.5")	10 per box
651	Medium, 13.7cm x 6.4cm (5.38" x 2.5")	10 per box
652	Small, 11.4cm x 4.4cm (4.5" x 1.75")	10 per box
653	X-Small, 10.8cm x 2.5cm (4.25" x 1")	10 per box

Dale[®]
www.dalemed.com

Dale[®] is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2017 Dale Medical Products, Inc.

PS-650-653 REV F