




# DNA-Free Cotton Tipped Applicator

Standard tip, Wood handle



## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free and RNase-free

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.188" Dia. (4.78 mm)
- 0.625" Length (15.88 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.62 mm)

## Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- Forensics / DNA testing
- Crime scene investigation
- Buccal cell collection
- Specimen collection



Ref #	Packaging
25-806 1WC FDNA Sterile	10/100/1 (10 boxes per case, 100 swabs per box, 1 swab per package)