

MANUMED EXERCISE COUCH



PARTNER FOR LIFE

Neurological exercise therapy with the Manumed Bobath. The large, extra wide treatment surface of 100 cm to 120 cm makes it easy to support and assist patients during therapeutic exercises. The Bobath and Vojta motion methods can be applied to perfection owing to the unmatched stability. A sturdy couch affording both patient and practitioner the most comfortable treatment positions.



- Very large surface (100 or 120x200 cm)
- Wheel raising mechanism for easy movement of the couch
- Electric or hydraulic height adjustment
- The 2-section model can be used as a neurological exercise or treatment couch



ENRAF-NONIUS B.V.

Vareseweg 127
3047 AT Rotterdam
The Netherlands
T: +31-(0)10 - 203 06 00
E: info@enraf-nonius.nl
www.enraf-nonius.com


MANUMED

MANUMED EXERCISE COUCH



PARTNER FOR LIFE


MANUMED EXERCISE 2-SECTION HEIGHT ADJUSTABLE



Type

- 511 Manumed Exercise 2-section, hydraulic H/L
- 521 Manumed Exercise 2-section, electric H/L


MANUMED EXERCISE 1-SECTION HEIGHT ADJUSTABLE



Type

- 510 Manumed Exercise 1-section, hydraulic H/L
- 520 Manumed Exercise 1-section, electric H/L

MANUMED EXERCISE 1-SECTION FIXED HEIGHT



Type

- 500 Manumed Exercise 1-section, fixed height of 55 cm

TECHNICAL SPECIFICATIONS

Size tabletop (lxw):	200 x 100 cm for models 500, 510 and 520 optional 120 cm width for 510 and 520 203 x 120 cm for models 511 and 521
Height adjustment:	44-97 cm
Lifting time (min.-max.):	approximately 18 sec.
Lifting capacity:	250 kg for Type 520 and 521 (electric) 200 kg for Type 510 and 511 (hydraulic) 135 kg for Type 500 (fixed height)
Power supply:	120/230VAC 50/60 Hz
Current consumption:	2.0 A max.

NOTE: ALL-AROUND OPERATION IS AN OPTION.



MANUMED

MANUMED COLOURS

A matching Enraf-Nonius couch for every interior. As in other practices, a pleasant ambience is important in physiotherapy practices. Clients must feel comfortable in the room where they are being treated. As each practice is furnished differently and has its own colour scheme, we regard it as important for the treatment couches to fit in with your look and feel. That is why we can manufacture our treatment couches in several frame and upholstery colours. Our imitation leather upholstery contains no phthalates and is flame retardant – a unique combination that makes this material suitable for medical use!




NATURAL



811
volcano black



850
walnut




851
sand



886
leaf



888
snow



LIFESTYLE



410
anthracit



406
petrol



407
cherry



408
terracotta




409
skin




EXECUTIVE



012
black




213
grey




199
lavender



204
sirius blue



169
cornflower
blue



170
baby blue



Frame colours



Ral 7035
grey



Ral 9010
white

OPTIONS⁽¹⁾ MANUMED EXERCISE COUCH


UPHOLSTERY OPTIONS

Section Width	Standard 4 cm thick	DeLuxe 6 cm thick	Extra hard 4 cm thick	Extra soft 4 cm thick
100 cm ⁽¹⁾	s	x	x	x
120 cm ⁽²⁾	s	x	x	x

s = standard x = not available


(1) = standard for 500, 510, 520

(2) = standard for 511, 521, optional for 510, 520



FOAM DENSITY:

Standard: 38 Kg/m³





HEIGHT ADJUSTMENT

- A. Hydraulic (510/511)
- B. Electric (520/521) Footswitch (standard)
- C. Electric (520/521) All-around (optional)



LIFTING CAPACITY

- A. Electric - max. 250 kg (standard 520/521)
- B. Hydraulic max 200 kg (510/511)



OPTIONAL FEATURES

- A. Safety switch



(1) = NOTE: Please order options together with the couch because of modification of the couch during production

FIXATION STRAP

3445474 Fixation strap with safety lock 350 cm



3445474



ACCESSORIES MANUMED EXERCISE COUCH

STOOLS (AVAILABLE IN ALL COLOURS SHOWN BELOW)

3xxx010	Stool, height adjustable 40-53 cm
3xxx011	Stool, mobile, height adjustable 59-85 cm
3xxx012	Stool, mobile, provided with footrest, height adjustable 59-85 cm
3xxx013	Stool, mobile, anatomic seat, height adjustable 43-56 cm
3xxx020	Stool, mobile, saddle seat, height adjustable 59-85 cm



3xxx010

3xxx011

3xxx012

3xxx013

3xxx020


ROLLS AND CUSHIONS (AVAILABLE IN ALL COLOURS SHOWN BELOW)

3xxx002	Halfroll, flat bottom, small, especially for the ankle	28x9x11 cm
3xxx003	Halfroll, flat bottom, large	56x16x18 cm
3xxx023	Massage roll	50x Ø15 cm
3xxx004	Neck cushion, hourglass type	30x15 cm
3xxx005	Massage cushion	30x28x8 cm
3xxx015	Wedge shaped cushion	50x30x15 cm
1550813	Plug for nosehole, triangular shaped (only available in black)	16,5 x 21 cm
3445297	Wedge shaped, Kaltenborn, rubber model (only available in black)	22,5 x 10,5 x 8,5
3446161	Support bag with sand 1,2 kg ⁽¹⁾	15x25 cm
3446162	Support bag with sand 2 kg ⁽¹⁾	18x27 cm
3446166	Support bag with sand 4 kg ⁽¹⁾	25x35 cm
3446167	Support bag with sand 7 kg ⁽¹⁾	30x45 cm

(1) = only available in grey



EXECUTIVE



NATURAL

LIFESTYLE

