



Applicator Stick

Wood



Description

This wood applicator is made of Northern Maine White birch which is naturally biodegradable with a splinter free finish.

Dimensions

5.812" Overall Length (147.625 mm)
0.083" Dia. (2.108 mm)

Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

Applications

- Multi-use applicator stick
- For applying or stirring medicines
- Can be used in any setting: Medical, Health & Beauty, Industrial, Tattoo, Veterinary



Puritan®

| Ref # | Packaging |
|----------------------|--|
| 807 Non-Sterile | 20/1000 (20 boxes per case, 1000 per box) |
| 25-807 2W Sterile | 10/100/2 Bulk (10 boxes per case, 100 individually wrapped packages per box, 2 per pkg) |