



Junior Tongue Depressor

Wood



- Biodegradable wood handle

Description

This junior size tongue blade is made of Northern Maine white birch which is naturally biodegradable. Available sterile or non-sterile depending on the application. Primarily used in oral patient care and as a first aid kit component.

Dimensions

5.5" Overall Length (139.7 mm)
0.625" Width (15.875 mm)

Qualities

- Medical grade quality
- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges
- Made in USA

Applications

- Oral patient care and examinations
- First aid kit component
- Veterinary examinations
- Applying or mixing medicines
- Health and beauty applications - waxing

Ref #	Packaging
711 Non-sterile	10/500 (10 boxes per case, 500 per box)
25-711 Sterile	10/100/1 (10 boxes of 100 individually wrapped-packages, 1 p/pkg)
712 Non-sterile	12/100 (12 boxes per case, 100 per box)



Puritan®